



LECKEY



**Say hello to the
Leckey family
of products.**

I am delighted to introduce our new Leckey product guide and celebrate our Leckey family circle.

At Leckey we believe in ‘Supporting Children’.

To everyone in Leckey, Supporting Children means working together with our family of experts – designers, therapists, parents and the broader clinical community to achieve our goals.

Our **Mission** is; To improve the quality of life of children with special needs.

Our **Vision** is; A world where each child is included and they can access their full potential.

Our **Courageous Goal** is; To prevent hip migration in children with CP by 2025.

In Leckey we see hip migration as the single biggest life event that our children may experience and has the greatest impact on their quality of life and their carers from that moment on.

Everything we do in Leckey, in partnership with our family of experts, is aimed at achieving these goals.

Our approach is simple – lead through innovation and clinical focus.

Through early intervention, childhood right up to adulthood, Leckey’s experienced team of designers and clinicians work together with our family of experts to develop exciting innovative products that make a real difference to the quality of life for the children and their families.

We aim to develop clinically and technically superior solutions, products and programs, that improve the child’s Mobility, Ability and Participation.

All of our designs undergo extensive clinical and usability trials with physiotherapists, occupational therapists, users and their families.

They are supported by in-depth clinical research studies where we work with leading universities and clinicians internationally to ensure that all Leckey products deliver on the clinical needs of the healthcare professionals and the social needs of the user.

We are excited to be part of the Sunrise Medical family and together we are committed to our Mission, Vision and Courageous Goals and Improving People’s Lives. We want to work with you in partnership to build our family - a community of experts with shared goals and values.

Join our family and together we will make a difference,



Ben Stocks
Managing Director



Contents

Seating and Positioning		04
Early Activity System		06
Pal Seat		08
Squiggles Seat		09
Mygo Seat		10
Everyday Activity Seat		12

Standing and Gait Trainers		14
Squiggles Stander		16
Squiggles+		18
Mygo Stander		20
Horizon		22
Pronestander		24
Freestander		25
Totstander		25
MyWay+		26

Other Products		28
Advanced Bath Chair		30
Shower Trolley		30
Sleepform		31



Seating and Positioning

LECKEY Early Activity System

Colorful, comfortable, supportive activity system for children aged 0-48 months.



Product Features

- Offers a gentle, playful introduction to 24-hour postural management.
- Practical all-in-one kit, providing versatile support. The positioning possibilities are endless and everything can be packaged up in the duffle bag for easy transport.
- Supplied with the popular Activity Development Programme which supports parent-led therapy at home.
- Soft-touch, anti-bacterial, machine washable and tumble dry safe covers.



Key Positions

- Back Lying [Supine]
- Tummy Lying [Prone]
- Side Lying
- Floor Sitting [Long sitting]
- Hands and knees [Four point kneeling]




LECKEY Pal Seat

Colorful, comfortable classroom seating system providing mild support for children aged 1-12 years.

Product Features

- Two-point lap belt and side pads for lateral stability. An optional ramped cushion can be added if extra help is needed to keep the pelvis and femurs in alignment.
- Where a less upright angle is required, the backrest angle and armrests adjust to suit the child's needs, and along with the large tray or table, provide appropriate support.
- Height adjustable armrests, legs and footplate.
- Available in two bright colors in wipe clean fabric (please see website for details).
- A range of base and footrest options available, including castors, rockers and stabilisers.



Technical Details				
	Size 1	Size 2	Size 3	Size 4
Age	1-3	2-5	4-8	7-12
User Height	Min 700mm [27.6"] Max 950mm [37.4"]	Min 900mm [35.4"] Max 1050mm [41.3"]	Min 1000mm [39.4"] Max 1300mm [51.2"]	Min 1200mm [47"] Max 1600mm [63"]
Max. User Weight	30 kg / 66 lbs	35kg / 77 lbs	40kg / 88 lbs	50kg / 110 lbs
Seat Depth	180mm-240mm [7"-9.4"]	220mm-290mm [8.6"-11.4"]	250mm-330mm [9.8"-12.9"]	300mm-400mm [11.8"-15.7"]
Seat Width	150mm-250mm [5.9"-9.8"]	170mm-270mm [6.6"-10.6"]	230mm-330mm [9"-12.9"]	280mm-380mm [11"-14.9"]
Seat Height	160mm-240mm [6.2"-9.4"]	185mm-285mm [7.2"-11.2"]	250mm-350mm [9.8"-13.7"]	300mm-400mm [11.8"-15.7"]

Supporting Children

LECKEY Squiggles Seat

Suitable for use as a home chair, school chair or with a mobility base, Squiggles meets all the postural and comfort requirements of a seating system for young children with special needs, while maintaining a fun, colorful design.

Product Features

- Designed specifically for the developmental needs of younger children- the one-piece seat base and adjustable sacral support can be used to accommodate flexion or encourage extension.
- The shape of the seat and back, the laterals with their fabric tip which wrap securely around the child and the way in which the 4-point pelvic harness integrates with the hip pads, provides a snug, comfortable fit and streamlined look.
- Being low to the ground, the Squiggles Seat promotes peer to peer interaction, making it ideal for home and nursery/school environments.
- Adjustable footplate can be tilted to the preferred angle, ensuring the feet are well supported and providing a secure base for maximum upper body function. The footplate flips away to allow for easy transfers.
- Tilt-in-space feature provides a simple way to achieve a change of position without disturbing overall posture.



Technical Details	
One Size Only	
Age	1-5
User Standing Height	750mm-1100mm [29.5"-43.3"]
Max. User Weight	22kg [48.4lbs]
Seat Depth	190mm-290mm [7.5"-11.4"]
Seat Width	160mm-260mm [6.3"-10.2"]
Chest Width	160mm-220mm [6.3"-8.6"]
Backrest Height	320mm-380mm [12.6"-15"]
Lower Leg Length	122mm-270mm [4.8"-10.6"]
Tray Height	115mm-170mm [4.5"-6.7"]
Backrest Angle	Prone 10° - Recline 25°
Seat Height (Hi-Low Chassis)	270mm-650mm [10.6"-25.6"]

LECKEY Mygo Seat

The Mygo offers postural support while enabling function and mobility. It is available in two sizes for children aged 3-14 years.



Product Features

- Pelvic Cradle with independent hip guides provides more proximal positioning for children with high tone. The flexible adjustment in the back section can encourage a posteriorly tilted pelvis into neutral alignment, and the overall design prevents users sliding forward.
- The Mygo Seating System's innovatively designed split seat base enables the legs to move independently beyond the midline to accommodate windsweeping postures. A range of up to 13° abduction and 10° adduction may be achieved whilst offsetting the hip guides can aid in accommodating more extreme windsweeping by providing additional space.
- The multi-adjustable 4-point pelvic harness on the Mygo Seating System ensures stability of the pelvis by providing support across the anterior superior iliac spines (ASIS). The 4-point attachments (two at the hip guides and two at 90°) mean that the harnesses stay in the correct anatomical position.



Available in 4 fabric colors

Technical Details		
	Size 1	Size 2
Age (approx)	3-10	8-14
User Weight	Min 18kg (40lbs) Max 50kg (110lbs)	Min 18kg (40lbs) Max 60kg (132lbs)
User Height	Min 1050mm (41.3") Max 1500mm (59")	Min 1270mm (50") Max 1680mm (66.1")
Seat Width	Min 200mm (8") Max 325mm (12.8")	Min 220mm (8.7") Max 345mm (13.6")
Seat Depth	Min 270mm (10.6") Max 420mm (16.5")	Min 350mm (13.8") Max 470mm (18.5")
Chest Width	Min 170mm (6.7") Max 270mm (10.6")	Min 170mm (6.7") Max 270mm (10.6")
Backrest Height	Min 360mm (14.2") Max 470mm (18.5")	Min 460mm (18.1") Max 570mm (22.4")
Backrest Angle	Prone 10° - Recline 25°	Prone 10° - Recline 25°
Lower Leg Length (Seat to Sandal)	Min 215mm (8.5") Max 350mm (13.8")	Min 315mm (12.4") Max 470mm (18.5")
Armrest Height	Min 160mm (6.3") Max 210mm (8.3")	Min 210mm (8.3") Max 260mm (10.2")
Tray Size	550 x 480mm (21.6" x 18.9")	550 x 480mm (21.6" x 18.9")
Seat Height (footpedal chassis)	340mm-655mm (13.4"-25.8")	340mm-655mm (13.4"-25.8")

LECKEY Everyday Activity Seat

Strong, durable, practical and adjustable seating system, providing moderate support for kids and young adults aged 1-18 years.



Product Features

- Ideal for everyday use- the soft-touch breathable fabric, well-padded ramped seat cushion, contoured back cushion and padded arm rests combine to offer the ultimate in comfort.
- The wide range of accessories and options make it easy to adjust the seat as the user grows and their needs change.
- The flip-up footplate, hi-low base (to floor), long armrests and 10° prone tilt-in-space support self-transfers.
- 4-point pelvic harness gives a secure, stable base of support.
- Height and angle adjustable chest laterals available.
- Easy care – covers can be removed and machine washed.



 Available in blue vinyl

Technical Details	Size 1		Size 2		Size 3	
	mm/kgs	inches/lbs	mm/kgs	inches/lbs	mm/kgs	inches/lbs
Max. User Weight	30kg	66lbs	50kg	110lbs	100kg	220lbs
Weight - Seat & Hi-Low Chassis	29.5kg	65lbs	34kg	75lbs	49kg	108lbs
Seat Height (Hi-Low Chassis)	250-600	9.8-23.6	320-600	12.6-23.6	450-600	17.7-23.6
Seat Width	160-280	6.3-11	260-380	10.2-15	340-500	13.4-19.7
Seat Depth (Standard)	150-280	5.9-11	230-380	9-15	330-480	13-18.9
Lateral Distance - Chest	150-220	5.9-8.7	170-240	6.7-9.4	240-370	9.4-14.6
Back Height	300-400	11.8-15.7	400-520	15.7-20.5	450-630	17.7-24.8
Back Angle	30° Supine 10° Prone		30° Supine 10° Prone		30° Supine 10° Prone	
Armrest Height	150-250	5.9-9.8	200-300	7.9-11.8	180-325	7-12.7
Lower Leg Length	150-250	5.9-9.8	250-450	9.8-17.7	250-450	9.8-17.7
Lower Leg Angle	-10° to +50°		-10° to +50°		-10° to +50°	
Footplate Angle	-15° to +10°		-15° to +10°		-15° to +10°	
Tilt-in-Space	30° Supine 10° Prone		30° Supine 10° Prone		30° Supine 10° Prone	



Standing and Gait Trainers

LECKEY Squiggles Stander

Colorful, comfortable and supportive standing system for children aged 1-5 years.



Product Features

- 3-in-1 prone, supine and upright stander.
- Ideal for home use due to ease of adjustments.
- Minimalist child-focused design with fun bright pads blends easily into nursery/school environments.
- Low to the ground, encouraging socialization and interaction at eye level.
- Individually adjustable chest and pelvic supports are cushioned for comfort and adjust easily in height, width and angle to suit individual needs, ensuring a comfortable and secure position.
- Frame can be easily disassembled for storage and transportation.



Available in 3 fabric colors

Technical Details	
One Size Only	
Age	1-5
User Height	730mm-1100mm [28.7"-43.3"]
Max. User Weight	22kg [48lbs]
Chest Support Height [footplate to armpit]	550mm-820mm [21.6"-32.2"]
Chest Width [distance between laterals]	160mm-230mm [6.3"-9"]
Chest Support Angle Adjustment	+/- 30° 50mm [1.9"] - Depth
Hip Support [height from footplate to midbuttocks]	300mm-530mm [11.8"-20.8"]
Hip Width [distance between hip guides]	160mm-230mm [6.3"-9"]
Distance between Midline of Knees	140mm-215mm [5.5"-8.4"]
Footplate Angle	Plantarflexion 10° Dorsiflexion 10°
Tray Height [footplate to top of tray]	555mm-825mm [21.8"-32.4"] - Prone 480mm-825mm [18.8"-32.4"] - Supine
Tray Angle Adjustment	30° - Prone 40° - Supine

Supporting Children

LECKEY Squiggles+

The Squiggles+ retains the best of the Squiggles stander but also offers heaps of added benefits so little ones can continue to enjoy standing and be part of the action.



Product Features

- Anatomically-correct abduction to stabilize the joint and provide optimum loading for early hip joint development. Squiggles+ offers 60° of hip abduction bilaterally in prone, upright or supine. It promotes proper alignment in the hips, knees and feet as pivot points are in line with the hips.
- Increased comfort and support with wrap-around knee supports. The comfortable medial knee supports distribute pressure and are especially suited to abducted standing. The padded wraps maintain the knee in neutral alignment.
- Providing greater confidence with transfers- the larger hip pads improve the carer's and child's confidence and feeling of security when transferring. The improved pivot chassis facilitates full horizontal loading.
- Activity tray with improved height and adjustability.
- Height-adjustable, foldable push handle.



Available in 3 fabric colors
and a vinyl version



Technical Details	
	One Size Only
Age	1-5
Max User Weight	22kg/48.4lbs
User Height	Min 730mm / 28.5" Max 1100mm / 43.3"
Chest Support	Angle adjustment +/- 30° Depth adjustment 50mm / 2 inches
Chest Support Height (footplate to armpit)	Min 550mm / 21.5 inches Max 820mm / 32 inches
Chest Width (distance between laterals)	Min 160mm / 6.3 inches Max 230mm / 9 inches
Hip Support (height from footplate to mid buttocks)	Min 300mm / 11.8 inches Max 530mm / 20.8 inches
Hip Width (distance between hip guides)	Min 160mm / 6.3 inches Max 230mm / 9 inches
Distance between Midline of Knees	Min 140mm / 5.5 inches Max 215mm / 8.4 inches
Footplate Angle	Plantarflexion 10° Dorsiflexion 10°

LECKEY Mygo Stander

Strong, practical and adjustable standing system for kids aged 4-14 years.



Product Features

- 3-in-1 prone, supine and upright stander.
- Designed to improve upper limb function. The chest pad design replicates the shape of the rib cage which allows the shoulder girdle to freely protract and retract in prone, while in supine, the shoulder support promotes shoulder protraction. Furthermore, Mygo offers a lower tray height than most standers.
- Accommodates hip and knee flexion contractures of up to 25 degrees. The movement of the footplate and knee supports are linked which allows the user to stand with hips and knees flexed. This results in a more upright posture combined with the ability to support users whose contractures previously made them unsuitable for standing therapy.
- Use of the pelvic derotation belt puts the user's pelvis in a central and neutral position with minimum effort.



Available in 3 fabric colors

Technical Details

	Size 1	Size 2
Age	4-10	8-14
User Height	1080mm-1450mm [42.5"-57"]	1350mm-1700mm [53.1"-66.9"]
Max. User Weight	50kg [110lbs]	60kg [132lbs]
Chest Support Height (footplate to armpit)	760mm-1080mm [29.9"-42.5"]	980mm-1300mm [38.5"-51.1"]
Chest Width (distance between laterals)	160mm-280mm [6.2"-11"]	200mm-320mm [7.8"-12.5"]
Hip Support Height (footplate to mid-buttocks)	460mm-730mm [18.1"-28.7"]	620mm-875mm [24.4"-34.4"]
Hip Width (distance between hip guides)	220mm-330mm [8.6"-12.9"]	280mm-390mm [11"-15.3"]
Footplate Angle	Plantarflexion 10° Dorsiflexion 10°	Plantarflexion 10° Dorsiflexion 10°
Max. Knee Contracture Accommodated	25°	25°
Tray Height Range (from top of hip pad)	0mm-150mm [0"-5.9"]	0mm-150mm [0"-5.9"]
Tray Angle Adjustment	45° - Prone 45° - Supine	45° - Prone 45° - Supine

Supporting Children

LECKEY Horizon

Robust, durable 3-in-1 stander
for users aged 4 years to adult.



Product Features

- Accommodation of leg length discrepancy, scoliosis, lumbar lordosis, kyphosis, knee flexion contractures and asymmetrical head position. The frame's growth potential and the wide range of accessories accommodate changing need and potential asymmetry as the user grows.
- Perfectly sized for total dependent transfers as well as supine standing. The large landing space allows for safe removal and the addition of lift slings coupled with support pads provide a proximal fit for the perfect combination.
- Robust, stable and durable frame ensures ongoing confident product use.
- Individually adjustable chest, hip and knee pads as well as footplates provide unsurpassed levels of adjustability.
- Height adjustable hip laterals and pelvic belt, enabling de-rotation of the hips, alignment of the pelvis and center of gravity over the base of support.



Available in grey

Technical Details		
	Size 2	Size 3
Age	9-18	14-adult
User Height	1400mm-1800mm [55.1"-70.8"]	1550mm-1800mm [61"-70.8"]
Max. User Weight	80kg [176lbs]	100kg [220lbs]
Shoulder Height (footplate to top of pad)	1000mm-1450mm [39.3"-57"]	1000mm-1450mm [39.3"-57"]
Hip Pad Height (footplate to top of pad)	750mm- 950mm [29.5"-37.4"]	750mm-1150mm [29.5"-45.2"]
Hip Pad Width	210mm-390mm [8.2"-15.3"]	310mm-470mm [12.2"-18.5"]
Knee Height (footplate to middle of pad)	400mm-550mm [15.7"-21.6"]	400mm-650mm [15.7"-25.5"]
Width Between Chest Laterals	240mm-380mm [9.4"-14.9"]	340mm-460mm [13.3"-18.1"]
Footplate Angle Adjustment +/-	10°	10°
Frame Angle Adjustment	0°-90°- Prone 0°-90°- Supine	0°-90°- Prone 0°-90°- Supine
Tray Angle Adjustment	0°-40°- Prone +/- 20°/30°- Supine	0°-40°- Prone +/- 20°/30°- Supine

Supporting Children

LECKEY Pronestander

Robust, adjustable prone standing system for users aged 1-18 years.

Product Features

- Flip-up sandals support users completing a standing transfer directly into the product.
- Long base ensures maximum stability for even the most active users, providing caregivers with peace of mind.
- Unique de-rotation belt and pelvic pad work together to keep the user in a central position within the frame with equal pressure over each foot.
- Individually adjustable chest, hip, knee and sandal supports to meet a wide range of clinical needs.
- Easy to set up and adjust with five incremental secure prone positions from 85 to 45 degrees.



Size 1 only available in Meadow Fabric
Size 2 & 3 only available in Mondeo Fabric

Technical Details			
	Size 1	Size 2	Size 3
Age	1-5	4-10	9-18
User Standing Height	800mm-1151mm [31.4"-45.3"]	1100mm-1448mm [43.3"-57"]	1400mm-1849mm [55.1"-72.7"]
Max. User Weight	60kg [132lbs]	80kg [176lbs]	90kg [198lbs]
Chest Pad Height	470mm-680mm [18.5"-26.7"]	680mm-975mm [26.7"-38.3"]	975mm-1350mm [38.3"-53.1"]
Prone Angles	85°, 75°, 65°, 55°, 45°		
Width between Chest Laterals	152mm-254mm [5.9"-10"]	203mm-305mm [7.9"-12"]	254mm-356mm [10"-14"]
Tray Sizes	540mm x 385mm [21.2" x 15.1"]	580mm x 410mm [22.8" x 16.1"]	645mm x 460mm [25.3" x 18.1"]

Supporting Children

LECKEY Freestander

Robust, upright standing system for users aged 1-18 years.

Product Features

- Flexibility to alter the frame to meet individual user needs e.g. pelvic support can be removed for users with trunk control, conversely, the tray height and angle can be increased where more support is required. A simple, active set-up is also possible for more able users.
- Straightforward, intuitive adjustment as all supports run off two vertical uprights.
- Pelvic belt achieves a unique position which enables de-rotation of the hips, alignment of the pelvis and center of gravity over the base of support as well as accommodation or correction of flexion at the hips and knees.



Size 1 only available in Meadow Fabric. Size 2 & 3 only available in Mondeo Fabric.

Technical Details			
	Size 1	Size 2	Size 3
Age	1-5	4-10	9-18
User Standing Height	800mm-1151mm [31.4"-45.3"]	1100mm-1448mm [43.3"-57"]	1400mm-1849mm [55.1"-72.7"]
Max. User Weight	60kg [132lbs]	80kg [176lbs]	90kg [198lbs]
Chest Pad Height	470mm-680mm [18.5"-26.7"]	680mm-975mm [26.7"-38.3"]	975mm-1350mm [38.3"-53.1"]
Width between Chest Laterals	152mm-254mm [5.9"-10"]	203mm-305mm [7.9"-12"]	254mm-356mm [10"-14"]
Tray Sizes	540mm x 385mm [21.2" x 15.1"]	580mm x 410mm [22.8" x 16.1"]	645mm x 460mm [25.3" x 18.1"]

LECKEY Totstander

Simple and easy to use, upright standing system for kids aged 1-5 years.

Product Features

- Can be used as a stander or by therapists for standing assessments. Promotes active standing, movement and weight shift.
- Compact and easily transported to different environments e.g. from therapy centres to the home.
- Provides a gentle introduction to rehab equipment via a simple, easy-to-use, child-friendly design.

Technical Details	
	One Size Only
Age	1-5
Chest Band Height	470mm- 640mm [18.5"-25.1"]
Base Depth	450mm [17.7"]
Base Width	450mm [17.7"]
Max. User Weight	22kg [48 lbs]



Available in 1 fabric color combination

Supporting Children

LECKEY MyWay+

MyWay+ is a unique open frame gait trainer designed to support children in an upright, hands-free position. This improves stepping to maximize opportunities for exploration, interaction, and participation.



Product Features

- The unique corset-like harness gives core trunk stability for optimal stepping pattern and active movement. The MyWay+ harness provides the support required for users to develop head and limb control. Its hug-like effect increases sensory feedback, dampens involuntary movement and increases confidence and the chance of succeeding.
- MyWay+'s open-front, hands-free frame maximizes opportunities for exploration and participation. The open-design also allows therapists to get close to the child and work with their lower limbs while they are in MyWay+, facilitating guided movement.
- MyWay+'s newly developed indoor/outdoor casters enable children to be active in lots of different environments. These casters are multifunctional and offer braking, direction lock, anti-roll-back protection, and resistance to assist steering, prevent backwards movement, and limit speed as required. The versatility of these casters allows the MyWay+ to be tailored to the individual child's needs, helping to maximize their opportunity for use and further promoting function. It is the only gait trainer that offers a side stepping feature.
- Range of optional accessories available to meet individual requirements, including adjustable headrest, narrow saddle, user handles, and ankle guides.



Narrow Saddle & Drop Attachment



3-Way Adjustable Headrest



3-Way Adjustable User Handles



Harness available in 6 accent colors

Technical Details			
	Size 1 Frame	Size 2 Frame	Size 3 Frame
Age [approx.]	1-5	4-11	9-16
User Shoulder Height [to floor]	590mm-870mm [23.2"-34.3"]	870mm-1190mm [34.3"-46.9"]	1070mm-1490mm [42.1"-58.7"]
User Saddle Height [inside leg]	240mm-410mm [9.4"-16.1"]	400mm-610mm [15.8"-24"]	550mm-810mm [21.7"-31.9"]
Max. User Weight	25kg [55.1 lbs]	50kg [110.2 lbs]	80kg [176.4 lbs]
Upper Support Unit Weight	Size 1 frame & upper support unit do not separate		11kg [24.2 lbs]
Frame Width	570 mm [22.4"]	670mm [26.4"]	720mm [28.3"]

Harness						
	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6
Age [approx.]	1-4	2-6	5-9	8-12	12-15	Up to 18
Waist at Navel Level	460mm-560mm [18"-22"]	500mm-600mm [19.5"-23.5"]	550mm-650mm [21.5"-25.5"]	600mm-700mm [23.5"-27.5"]	700mm-800mm [27.5"-31.5"]	800mm-900mm [31.5"-35.5"]
Color	Purple	Green	Blue	Red	Orange	Grey



Other
Products

LECKEY Advanced Bath Chair

The Leckey Advanced Bath Chair provides postural support while ensuring the individual is safe and secure during bathing.

Product Features

- Provides the user with appropriate levels of postural support ensuring they are safe and secure at all times during bathing.
- Lightweight and robust frame.
- Simple and easy-to-use adjustments mean that the product can be set up in a variety of positions in seconds.
- Height adjustable legs.
- Breathable cover, removable and machine washable.
- Back angle adjustment from 0-90 degrees.
- Adjustable hip belt.
- Chest belt with integrated lateral supports.
- Easy to use hand adjustments.
- Folds flat for storage.



Available in 2 fabric colors

Technical Details				
	Size 1	Size 2	Size 3	Size 4
Age	1-5	4-9	8-14	12-18
Min User Height	750mm [30"]	950mm [37"]	1150mm [45"]	1350mm [53"]
Max User Height	1050mm [41"]	1250mm [49"]	1450mm [57"]	1650mm [65"]
Max User Weight	72kg [159lbs]	72kg [159lbs]	72kg [159lbs]	72kg [159lbs]

LECKEY Shower Trolley

The Shower Trolley complements the Advanced Bath Chair which can be converted within seconds into a shower chair.

Product Features

- Converts the Advance Bath Chair to a versatile shower chair in seconds.
- Manufactured from medical grade stainless steel.
- Robust frame and high quality lockable castors.
- Can be easily and safely manoeuvred over wet and slippery surfaces.
- Can be adjusted to fit all four sizes of the Leckey Advance Bath Chair.



Technical Details	
	One Size Only
Age	1-18
Length	900mm-1150mm [36"-45"]
Width	600mm [24"]
Seat Height from Floor	650mm-700mm [26"-28"]

LECKEY Sleepform

Strong, durable, practical and adjustable sleeping system for kids aged 0-18 years.

Product Features

- Close proximal fit cocoons the child and makes them feel safe and secure. Adapts quickly and easily to changing body shape.
- Compact and portable for visits to family, going on holiday etc. Suitable for use on a variety of surfaces.
- Supports 24-hour postural management as can be used during the day to support changes of position.
- Can be used in cots and single beds and suits all sleeping positions.
- Air flow mattress helps with temperature control and comfort.
- Grid sheet to record best sleeping positioning.
- Temperature control sheet available.



Technical Details				
	Size 1	Size 2	Size 3	Size 4
Age	0-1	1-5	1-5	5-18
Suitable for	Crib	Crib	Single Bed	Single Bed

Supporting Children



Sunrise Medical (US) LLC, Fresno, CA 93727
MK-129916 Rev. B ©12.2023

800.333.4000 www.SunriseMedical.com

