



Syncra™ Seating

Exclusively for the QUICKIE® SEDEO® PRO™

A cohesive and attractive seating system specifically designed for tilt and recline.



JAY Syncra

Exclusively for the QUICKIE® SEDEO® PRO™.

JAY® Syncra™ is a cohesive and attractive seating system specifically designed for optimal performance with tilt and recline systems. Syncra has a width and height adjustable back, a cushion that fits mild to complex seating needs, and positioning accessories for support and stability.

JAY® Syncra™ is the optimal seating for power wheelchair users who have a history of pressure injuries.

JAY Fluid with Cryo™ Technology



The industry favorite JAY® Fluid is now enhanced with the revolutionary, patented CRYO™ TECHNOLOGY.

The CryoFluid™ insert addresses the four main risk factors for tissue breakdown: pressure, shear, temperature, and moisture.

CryoFluid™ actively cools the skin for up to eight hours*, lowers skin temperature and reduces sweating, therefore reducing moisture concerns.

JAY® Fluid evenly distributes pressure across the seating surface and around bony prominences.

The CryoFluid™ PLA insert is divided into five chambers, the two chambers under the ischial tuberosities (ITs) include CryoFluid and the other three chambers include industry favorite, JAY Flow™ Fluid.

Learn more at: www.SunriseMedical.com/Cryo

*Internal testing data at 25°C. Results may vary.



Positioned for Recline

Syncra™ Cushion

- Engineered to accommodate pelvic movement through recline
- Lowered rear pelvic wall
- Deeper leg troughs

Syncra™ Backrest

- Scalloped back shell allows better posture during recline
- Optimal shell placement ideal for recline positioning
- Built-in adjustability – including 2" of width and height – accommodates patient changes

Positioners & Supports

- Mounting positions optimized for WHITMYER® headrests with 2" lateral movement and left, right, or center mounting options
- Designed for integrating multiple lateral, trunk, and thigh supports



Thigh Supports



The Syncra™ removable lateral thigh supports are height and angle adjustable to align with the angle of the leg for optimal positioning. The quick-release lever allows for easy placement and removal, even for clients with limited dexterity.

Syncra Cushion Options

The Syncra™ cushion has multiple covers and inserts that can be interchanged based on the client's needs.

Cushion Inserts



JAY Fluid with Cryo Technology



Visco Foam



Air Insert

Cover Options



Premium Leatherette



4-Way Stretch



3DX Microclimatic