

AI
ALUMINIUM 7020



CE

Customer Contact Details

Customer Name: _____

Advisor: _____

Address: _____

Phone: _____

Email: _____

Ship to Customer Address:

Shipping Name: _____

Shipping Address: _____

Phone: _____

Email: _____

Additional Information

Notes: _____

Current Wheelchair: _____

Customer Weight: _____

Customer Height: _____

Customer Disability: _____

Other

VAT Exemption Declaration (UK Only):
I declare that I am chronically sick or disabled by reason of:

and that I am receiving from RGK Wheelchairs Ltd, the wheelchair specified herein and I claim relief from Value Added Tax under group 12 of schedule 8 to Value Added Tax Act 1994.
Signature: _____ Date: _____

RGK Contact Authorisation:
From time to time, RGK would like to send you exclusive information, news and offers on everything RGK, by email or post. We will not pass any of your details to any third party.
I would like to be contacted by RGK:
Signature: _____ Date: _____

SERIAL NO: _____

Guide:

- Base product
- Required option
- Optional Extra
- Critical Information
- Information Needed

PRODUCT STYLE

EXAXZ010 ELITE X
\$5,847.00

XA600 Heat treated frame

EXAXA062 Standard

EXAXA063 Standard with 25mm seatrail
\$69.00



Max User Weight 125kg (275lbs)

PLEASE NOTE: ALL RGK CHAIRS ARE HAND BUILT AND HAVE A MANUFACTURING MEASUREMENT TOLERANCE OF +/- 5MM

CAD APPROVAL

XW009 QUOTE ONLY (NO CAD FOR APPROVAL)
\$0.00

XW011 ORDER (NO CAD FOR APPROVAL)
\$0.00

XW012 QUOTE CAD drawing for approval
\$421.00 (2 customer amendments ONLY)

XW010 ORDER CAD drawing for approval
\$260.00 (2 customer amendments ONLY)

KEY MEASUREMENTS

XY010 Seat width (Inside to inside of sideguards)

(Max - 500mm)

XY020 Seat length

(Max - 500mm)

XY030 Extra seat length (Provides additional leg support & hand grip)

XY050 Front frame height

(Max - 630mm)
(Min - 360mm)

XY060 Rear frame height


(Max - 630mm)
(Min - 360mm)

XY070 Frame length

(Min = XY020 + XY030 + 140mm)

XA020 Ergonomic seat
\$80.00

Specify length (Standard = 120mm)



FRONT FRAME AND SEAT

Front Frame Options

XY080 Seat taper design


i Specify taper width mm



XA130 V-seat design
\$212.00

i Specify taper width mm


NOTE: INSIDE TO INSIDE WHEEL (XY470) REQUIRED



Advanced Front Frame Options

XY082 Seat taper position
\$83.00

i Specify measurement mm




XA131 V-seat and taper design
\$286.00

i Specify inside width at front of sideguard (A) mm


i Specify taper width (B) mm

NOTE: INSIDE TO INSIDE WHEEL (XY470) REQUIRED



XA136 Parallel taper design
\$83.00 (Not with Seat Taper XY080)

i Specify taper width mm




XA135 Asymmetrical taper design
\$351.00

i Specify taper width (A) mm

i Specify offset (B) (Inside of sideguard to inside of taper width) mm

Specify offset taper side

i Left Right




Seat Options

XR050 Carbon fibre seat plate
\$337.00 (riveted to frame)

i Specify length of plate mm (Std XY020 - 20mm)

XA260 Drop seat design
\$194.00 (only with either XR040 or 050)

i Specify drop mm (Std = 20mm)



FOOTREST

Footrest Type


EXAXB080 Fixed position footrest

EXAXB145 Multi position footrest tubes
\$140.00



Footrest Measurements

Front footrest height
[] mm (min = 35mm)



Footrest width
[] mm (Std = 200mm)





Rear footrest height
[] mm (min = 35mm)



Footrest length
[] mm (Std = 250mm)

Footrest position
(front of footplate to castor pot)
[] mm



Footrest Plate

XB190 Carbon fibre footplate
\$219.00

XB100 ABS moulded foot tray
\$129.00

EXAXB140 Spacer kit
\$19.00 (Allows for 25mm increase in height)

Specify length
[] mm (Max = 380mm (Full)
Std = 330mm)

Footrest Plate Options

XB240 Aluminium raised edges
\$75.00 (Cost per Edge)

Front Back Left Right

Specify height
[] mm (Max = 60mm)

XB250 Carbon Fibre raised edge
\$290.00 (1 side only)

Front Back


Specify height
[] mm (Max = 45mm)

Castor Pin Width

XY330 Castor pin width

[] mm **MAX 320mm**

Optimise measurements



CAMBER BAR, WING AND ANTI-TIP

Centre of Gravity and Camber Angle

EXAXY510 Fixed centre of gravity

Specify measurement
[] mm

(Note: dimension is taken from front of backrest to centre of axle)



XY390 Camber angle

[] 18 [] 20



Wing Options

EXAXA110 Round design basketball wing

EXAXA125 Arrow design basketball wing
\$176.00

EXAXA120 Straight sided design basketball wing
\$176.00

EXAXA129 Individual design basketball wing
\$176.00

Anti-tip options

EXAXJ040 Single fixed anti tip

EXAXJ060 Double fixed anti tip
\$249.00

SIDEGUARDS

Sideguard Style


EXAXE230 Welded aluminium sport
\$374.00

XE310 High design upgrade
\$176.00




EXAXE250 Welded aluminium sport with Carbon Fibre
\$738.00


Sideguard Design (Choose One)



Straight design




Round design



Follow tyre

Sideguard Measurements

Sideguard length

 mm


Sideguard height (back)

 mm


Sideguard height (front)
(optional)

 mm



Follow Tyre Design

25mm above tyre

Individual above tyre

 mm

Level with tyre




BACKREST

Backrest Height

XY230 Fixed backrest height



Specify height


 mm


Ergonomic Options

XY360 Backrest angle

-5 -3 0 +3 +5

Individual Deg
(Max -10/+10 Deg)





XG100 Ergonomic backrest
\$80.00



Vertical height mm




Rearward position mm



XG110 Inward taper
\$80.00
(inside of backrest posts)




Specify width / per side mm



XG120 Outward flare
\$80.00
(inside of backrest posts)



Specify width / per side mm



CASTOR AND REAR WHEELS

Castor Options

- | | |
|--|---|
| <input type="checkbox"/> XC160 Sport 61mm castor (Black) | <input type="checkbox"/> XC181 Bones Reds bearings (All castor bearings) |
| <input type="checkbox"/> XC200 Sport micro fork | <input type="checkbox"/> XC182 Bones Swiss Ceramic bearings (All castor bearings) \$278.00 |

Rear Wheel Size

- | | | | |
|---|---|---|--|
| <input type="checkbox"/> XY420 24" Wheels | <input type="checkbox"/> XY430 25" Wheels | <input type="checkbox"/> XY440 26" Wheels | <input type="checkbox"/> XY450 700C Wheels |
|---|---|---|--|

Rear Wheel Gap

<input type="checkbox"/> XY460 Rear wheel gap <input type="checkbox"/> 15 <input type="checkbox"/> 20 <input type="checkbox"/> 25 <input type="checkbox"/> 30 Individual <input type="text" value="mm"/>	OR	<input type="checkbox"/> XY470 Inside to inside of wheel <input type="text" value="mm"/> NOTE: REQUIRED FOR V SEAT AND TAPER
--	----	--

Rear Wheel Options

- | | |
|--|---|
| <input type="checkbox"/> XK080 Spinergy SLX 24 spoke \$750.00 | <input type="checkbox"/> XK330 Spinergy SLX R10 heavy duty 24 spoke (18 deg camber angle only) \$1,385.00 |
| <input type="checkbox"/> XK101 Spinergy XSLX 24 spoke + Caps \$822.00 | <input type="checkbox"/> XK521 Spinergy Ceramic bearing upgrade \$252.00 |
- Spinergy Colour Options:**
- | | | | |
|--|--|--|---|
| <input type="checkbox"/> Black spokes | <input type="checkbox"/> White spokes | <input type="checkbox"/> Yellow spokes | <input type="checkbox"/> XK180 Individual spoke design (Specify colour) \$169.00 |
| <input type="checkbox"/> Red spokes | <input type="checkbox"/> Blue spokes | <input type="checkbox"/> Pink spokes | Inside: _____ |
| <input type="checkbox"/> Purple spokes | <input type="checkbox"/> Orange spokes | <input type="checkbox"/> Green spokes | Outside: _____ |
- Spinergy Anodise Options**
- | | | | |
|--|--|--|---|
| <input type="checkbox"/> XK190 Black hub (std option) | <input type="checkbox"/> XK200 Silver hub \$92.00 | <input type="checkbox"/> XK201 Red hub/rims \$274.00 | <input type="checkbox"/> XK205 Individual hub/rims Select colour: _____ \$316.00 |
| <input type="checkbox"/> XK202 Purple hub/rims \$274.00 | <input type="checkbox"/> XK203 Blue hub/rims \$274.00 | <input type="checkbox"/> XK204 Gold hub/rims \$274.00 | |

Supply Without Wheels

- | | | |
|---|----------------------------------|---------------------------------|
| <input type="checkbox"/> XK040 Supply without Wheels -\$184.00 | Wheel Type: <input type="text"/> | Tyre Type: <input type="text"/> |
|---|----------------------------------|---------------------------------|

Handrim Position / Options

- | | |
|---|--|
| <input type="checkbox"/> XY480 Narrow set handrim | <input type="checkbox"/> XY485 Wide set handrim |
| <input type="checkbox"/> XL090 Remove excess mounting tabs (recommended for small wheel gaps) \$101.00 | <input type="checkbox"/> XL080 Super narrow set handrim (only with XL040 and includes XL090) \$203.00 |

Handrim Style

- | | |
|---|---|
| <input type="checkbox"/> XL010 Aluminium hard anodised 6 tab | <input type="checkbox"/> XL040 Titanium 6 tab (Not with 700C) \$240.00 |
| <input type="checkbox"/> XL070 Aluminium hard anodised 12 tab (24" not available) \$101.00 | <input type="checkbox"/> XL050 Stainless steel 6 tab \$132.00 |

Tyre Options

- | | |
|--|--|
| <input type="checkbox"/> XM070 Kenda Kaliente <input type="checkbox"/> Red <input type="checkbox"/> Blue | <input type="checkbox"/> XM095 Kenda Wide (Grey Only) \$35.00 |
| <input type="checkbox"/> XM150 IRC Exeracer Pro (Red Only) \$83.00 | |
| <input type="checkbox"/> XM080 Tufo MS3 tubeless (Blue only) \$83.00 | |

Spare Wheel

- | | |
|--|--|
| <input type="checkbox"/> EXAXK500 Spare wheel Spinergy (Built up + axle pin) \$845.00 | <input type="checkbox"/> XK523 Spare wheel upcharge - Individual Hubs/Rims \$122.00 |
| <input type="checkbox"/> XK520 Spare wheel upcharge - Titanium Rims \$196.00 | <input type="checkbox"/> XK524 Spare wheel upcharge - Individual Spoke Colour \$88.00 |
| <input type="checkbox"/> XK522 Spare wheel upcharge - Remove Excess Tabs \$51.00 | |

Inner Tube Valve Options

- | | |
|---|---|
| <input type="checkbox"/> XY500 Schraeder inner tube valve | <input type="checkbox"/> XY490 Presta inner tube valve (Standard if no option selected) |
|---|---|

FINISH OPTIONS

Frame finish

XN010 Powder coated frame

Frame Comp

Gloss

Frame Comp

Matt

Frame

Special

- Black gloss 103
- White gloss 100
- Grey gloss 132
- Cream gloss 135
- Red gloss 108
- Yellow gloss 120
- Forest green gloss 143
- Green gloss 107
- Dark blue gloss 151
- Royal blue gloss 154
- Blue gloss 174
- Duckegg blue gloss 141
- Pink gloss 188
- Orange gloss 180
- Purple gloss 157

- Black matt 105
- White matt 101
- Grey matt 131
- Light grey matt 133
- Navy blue matt 152

- Chrome Mirror 300
- Sparkle silver 170
- Metallic black matt 104
- Metallic grey matt 130
- Black sandpaper 201
- Grey sandpaper 202
- White sandpaper 200

XN200 Individual paint colour (State Colour/RAL code)
\$172.00

EXAXN030 Brushed frame
\$347.00

XN160 Brushed components
\$131.00
(not required with fixed footrest & backrest)

EXAXN050 Polished frame
\$919.00

XN170 Polished components
\$193.00
(not required with fixed footrest & backrest)

Accent Colour

Note: Accent colour includes, forks, castor caps, decals, ELITE logo and stitching. Cushion, toe strap, frame protectors and RGK logo are

- Black accent
- Red accent
- Purple accent

- Silver accent
- Blue accent
- Carbon accent
(caps only, black for others)

XN250 Individual accent
\$172.00

Select Colour:

UPHOLSTERY

Seat Upholstery

EXAXP040 Seat sling strap adjustable

No seat sling (Seat Plate)

Backrest Upholstery

EXAXP055 Airtech sport backrest strap adjustable
(comes with Ti buckles)

XP070 Backrest upgrade to leather corners
\$120.00
(Black)

Backrest Stitching Options

- Silver
- Black
- Blue
- Individual:
- Red
- Orange
- Purple

RGK Logo Colour

- White (Std)
- Green
- Blue
- Individual:
- Red
- Orange
- Purple

Model Logo Colour

- White (Std)
- Green
- Blue
- Individual:
- Red
- Orange
- Purple

Backrest Accessories

XP170 Individual backrest embroidery
\$136.00
ADD DETAILS ON PAGE 9

CUSHION

Cushion foam options

- XP080** \$83.00 Cushion foam
(select density and size below)
 - Soft (CMS040) Medium (CMS050)
 - 50mm 75mm 100mm Other mm
- XP340** \$78.00 Hard Plastazote cushion upgrade
(Indicate max height)
 - 580mm 630mm

Cushion Cover Upgrades

- XP085** Nylon cushion cover (Black)
- XP090** \$64.00 Air Lite cushion cover (Black)
- XP380** \$89.00 Neoprene cushion cover (Black)

STRAPS

Velcro straps

- XS050** \$119.00 100mm Pro adjustable velcro strap
- XS055** \$119.00 50mm Pro adjustable velcro strap
- XS057** \$119.00 50mm Pro adjustable velcro **Knee** strap
(measurement from tube to tube over knees)

Velcro strap position

Strap 1 Type

Specify strap length mm

Specify distance from back post mm

Strap 2 Type

Specify strap length mm

Specify distance from back post mm

Ratchet Straps

- XS020** \$81.00 Sport ratchet strap (Click strap)
(11" x 2" pad)
- XS021** \$81.00 Sport ratchet strap (Single)
- XS025** \$212.00 Pelvi Loc 2DA ratchet strap
(For thighs and knees)
Select size:
 - Medium (12") Large (14")
- XS030** \$227.00 Pelvi Loc 3DA ratchet strap
(For hips)
Select size:
 - S (8") M (12") L (14") XL (17")
- XS035** \$135.00 M2 ratchet strap
(For thighs and knees)
 - 8" 10" 12"
- XS036** \$135.00 M2 ratchet strap (Single)
(For Feet)
 - 8" 10" 12"

Ratchet strap position

Strap 1 Type

Specify distance from back post mm

Strap 2 Type

Specify distance from back post mm

XS040 \$85.00 Fixed mounting points
(one set per click strap) x1 x2 x3

FRAME OPTIONS

- XR060** \$101.00 Calf strap mounting loops
- XR065** \$101.00 Strap location loops

ACCESSORIES

- XR030** \$99.00 Triple wheel bag (Black)
- XP130** \$90.00 Frame protectors full length
- XA150** \$64.00 Plastic tube fitted to offensive wing
- XS070** Foot strap *(select if required)*
- XS060** Calf strap *(select if required)*
- XP400** Back bar protector
- XR090** Tool kit

ADDITIONAL NOTES

INDIVIDUAL BACKREST EMBROIDERY

PLEASE PROVIDE JPEG IMAGE AT TIME OF ORDER AND NOTE POSITION BELOW



BRACE BAR

SPEC NOTES:

OFF SPEC REQUEST*

(please specify below the page number(s) and details to which the off spec request applies.

OFF SPEC REQUEST APPROVED (for internal RGK office only)

x

RGK Internal Sign-off

Off spec requests require a feasibility study which may increase order price and lead time. Any dimension/ option/ frame request not listed on the order form will be treated as an off spec request.*